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Equipment Safety Instruction: Hot air oven

Task	Hazard	Risk	Control
Switch on equipment	Electricity	- Workers electrocute - Spill liquid cause an electric shock	- Wear the proper PPE
Operation	Fire hazard	- Fire - Melting plastic can cause fires	- Do Not leave the oven running over night (If leaving oven running overnight/leave note for other lab employees listing: name and phone number, sample, time and temperature - In the case that the plastic is dry-heat sterilized in a hot air oven, the temperature should not over 160 °C - Ensure that the materials placed inside for heating are heat resistant material - Do not rinse glassware with organic solvents before being dried in hot air oven - Keep flammable and combustible materials away from hot air oven
Entry or exit of materials	Burn hazard	- Burns from heated glassware or touch the hot part of equipment	- Equipment operator training - Wear the proper PPE

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Equipment Safety Instruction : Hot air oven

Task	Hazard	Risk	Control
Entry or exit of materials	Glassware	- Workers injured from damaged or broken glass	- Do not bring broken glassware dried in hot air oven - Ensure the equipment is operating properly before bringing materials out - Do not overfill glassware or material
Maintenance			- Preparing a cleaning schedule

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