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Equipment Safety Instruction: Autoclave

Task	Hazard	Risk	Control
Before using the autoclave	Water level is not appropriated	- Damaged machine	- Check that water level in the autoclave is sufficient - Equipment operator training
	Temperature reset is not working	- Failed to properly sterilize equipment	- Routine inspection - Routine maintenance
Switch on equipment	Electricity	- Workers electrocute	- Wear the proper PPE
Entry or exit of materials	Working posture	- Workers injured	- Working position adjustment - Lifting an appropriate amount of weight
	Dropped materials	- Liquid spillage	- Spill kit - Remove items carefully
	Shattered glassware	- Workers injured	- Check the air pressure scale set at zero before opening
	Burn hazard	- Worker touch hot objects or touch hot parts of the machine	- Check the air pressure scale set at zero before opening - Open the door cautiously. Stand behind the door and slowly open it. Allow all steam to escape before reaching inside. - Wear the proper PPE
Autoclaved infectious materials	Hazardous biological agents	- Contaminate - Workers infected	- Wear the proper PPE - Do not overload the autoclave
Maintenance			- Preparing a cleaning schedule

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Equipment Safety Instruction: Autoclave

Practices for autoclaving liquids

1. Use only borosilicate glass (Pyrex™ or Kimax™) which can withstand the high autoclave temperatures.
2. Bottles of liquid should not be more than 2/3 full. To prevent bottles from shattering during pressurization and heating, the caps of containers with liquids must be loosened before loading.
3. Liquids should be within a heat resistant plastic tray (polypropylene or polycarbonate) containing an inch of water, to ensure even heating.

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