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Equipment Safety Instruction: Fume hood

Task	Hazard	Risk	Control
Before using the fume hood	Fume hood is not operating properly	<ul style="list-style-type: none"> - Generation of heat - Fire - Explosion 	<ul style="list-style-type: none"> - Smoke testing (Airflow velocity is not less than 100 ft/min) - Face velocity is not less than 0.4 m/min
During the operation	Chemical spills	<ul style="list-style-type: none"> - Contaminate - Skin burns or irritation caused by contact with the hazardous chemicals - Cancer occurring years after exposure to a carcinogenic substance - Chemical splash into an eye - Poisoning by absorption through the skin of a toxic substance - Accidental inhalation or ingestion 	<ul style="list-style-type: none"> - Review the SDS - Get everything out of the fume hood except what you need for your work - Work with the fume hood sash in the lowest possible position no higher than 18 inches - Keep chemicals and other materials at least 6 inches behind the plane of the sash - Don't block baffles or exhaust slots while working in the fume hood, this would prevent air from flowing into and out of the hood - Keep all containers of chemicals closed when not in use - Never remove the airfoil or modify chemical fume hood in any way - Keep all doors and windows closed while is going on in the fume hood - Move slowly in the hood - Wear the proper PPE - Preparation of a first aid kit, spill kit, eye wash and emergency shower

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Task	Hazard	Risk	Control
During the operation	Use of the incorrect hood or incorrect use of a hood	- Generation of heat - Fire - Explosion	- Do not use fume hood for long-term chemical storage - Never put your head into the hood while there is work in it
After the operation			- At the end of your experiment remove equipment and clean the surfaces - Dispose of chemical waste in a safe appropriate manner as identified by SDS and in accordance with laboratory rules
Maintenance			- Preparing a cleaning schedule

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