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Equipment Safety Instruction: Refrigerator 2-8 °C and freezer

Task	Hazard	Risk	Control
Before using the refrigerator	Storage of food and drink	- Contaminate	- Never store food in a refrigerator used to store chemical
	Improper installation	- Refrigerator will not be able to maintain the desired temperature	- Install the unit in a level area free from vibration with a minimum of 6 inches of space on the sides and rear and 12 inches at the top
	Electricity	- Workers electrocute	- Plug the unit into the emergency power back-up system and/or UPS - Wear the proper PPE
During the operation	Chemical spills Biological spills	- Chemical reaction : heat, explosion - Workers infected - Contaminate - Chemical splash into an eye - Breathe in hazardous chemical vapors	- Review the SDS and PSDS - All items stored are to be appropriately labeled - Wear the proper PPE
	Containers are not appropriate	- Liquid spillage	- Containers placed in a refrigerator or freezer should be completely sealed or capped - All liquid chemicals should be stored in secondary container - Avoid capping materials with aluminum foil, corks, and glass stoppers
	The door not close tightly	- Refrigerator will not be able to maintain the desired temperature	- Close the refrigerator immediately after putting in any items to maintain the internal temperature

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Equipment Safety Instruction: Refrigerator 2-8 °C and freezer

Task	Hazard	Risk	Control
During the operation	Storing flammable solvents Storing incompatible chemicals	Explosion	<ul style="list-style-type: none"> - Chemicals stored inside the refrigerator should be compatible (Flammable liquids may never be stored with strong oxidizers or acids) - Wrap caps of volatile flammable materials in parafilm and place its in sealed plastic bags
Maintenance			<ul style="list-style-type: none"> - Refrigerator should be defrosted at least once a year or as needed - Refrigerator should be cleaned using mild detergent and water at least once a year or as needed - The temperature will be recorded at least once daily - Refrigerators should be periodically inspected at least monthly

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