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Equipment Safety Instruction: Hotplate stirrer

Task	Hazard	Risk	Control
Switch on equipment	Electricity	- Workers electrocute	- Use a properly grounded electrical outlet with the proper voltage and current handling capacity
		- Spill liquid cause an electric shock	- Place Hotplate/Stirrer on flat and level surface - Wear the proper PPE
		- Fire and explosion	- Do not store combustible materials or flammable liquids near hot plates - Ensure that electrical cords do not come in contact with the hot plate surface
During the operation	Thermal	- Worker touch hot objects or touch hot parts of the machine	- Wear the proper PPE - Request stir and heat labels for hot plates when in use
	Electricity	- Fire and explosion - Spill liquid cause an electric shock	- Do not leave a hotplate stirrer unattended when in use - Wear the proper PPE
	Equipment turn on when not in use	- The liquid splashed into face and eyes	- Adjust speed from slow to fast to ensure sample is not stirred to quickly which can result in spills - Wear the proper PPE
	Chemical spills	- Working area is damaged - Workers injured	- In the event that volatile flammable and toxic materials are heated, there must be sufficient airflow in the fume hood to prevent a flammable/explosion - Provide secondary containment for any flammable liquids being heated to prevent liquid from contacting the hot plate in the event of a leak or overflow - Preparation of a spill kit

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2563

Task	Hazard	Risk	Control
During the operation	Glassware	- Workers injured	- Make sure the glass being heated is made of heat resistant material
Maintenance		- Contaminate	- Preparing a cleaning schedule - Unplugging when not in use - Regularly check hot plates for damage such as damaged electrical cords

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