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Equipment Safety Instruction: Vortex

Task	Hazard	Risk	Control
Switch on equipment	Electricity	- Workers electrocute	- Use a properly grounded electrical outlet with the proper voltage and current handling capacity
		- Liquid spill cause an electric shock	- Place Hotplate/Stirrer on flat and level surface - Wear the proper PPE
		- Fire and explosion	- Do not store volatile flammable materials near a centrifuge
During the operation	Electricity	- Fire and explosion	- Chemical/toxic substances must be vortexed inside a fume hood - Must not be used with explosive or highly reactive substances
	Container is not completely close	- Liquid spillage - Damage to vortex rotor	- Wear the proper PPE - Never overfill centrifuge tubes (do not exceed 3/4 full) - Preparation of a spill kit, first aid kit - Wear the proper PPE
	Work with hazardous substances	- Injury through exposure to aerosols with chemical and biological hazards	- Place a biohazard label on the vortex - Infectious and toxic chemical substances must be vortexed inside a biological safety cabinet and fume hood, respectively.
	-Chipped or cracked glassware	- Injury through exposure to aerosols with chemical and biological hazards	- Centrifugation at too high a speed for long time - Inspect all glassware to be used in the centrifuge (chipped/cracked) - Preparation of a spill kit, first aid kit
	Vortex at too high a speed for long time	- Liquid spill cause fire and explosion	- Slowly rotate the speed knob clockwise until the desired mixing level is achieved

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Task	Hazard	Risk	Control
Maintenance		- Contaminate	- Preparing a cleaning schedule - Never clean rotors and associated parts with abrasive wire brushes - Decontaminate the outside of the cups/buckets and rotors before and after centrifugation - Unplugging when not in use

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