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Equipment Safety Instruction: Biological Safety Cabinet

Task	Hazard	Risk	Control
Before using the biological safety cabinet	Contaminate	- Flow velocity is not appropriate	- Allow 5 minutes of operation to purge the system (check the flow velocity) - Clean the surfaces of the BSC with 70% alcohol before work
Entry or exit of materials	Hazardous biological agents spills	- Workers infected	- Wear the proper PPE - Preparation of a spill kit
	Chemical spills	- Irritation of eyes/skin with disinfectant chemical	- Do not use volatile chemicals - Preparation of a first aid kit
	Overfilling the cabinet with work materials Place any objects over the front air intake grill and block the rear exhaust grill	- Disturbance of inflow air velocity - Contaminate	- Do not obstruct or block any of the air grilles in the front or back of the work zone. - Keep only necessary equipment or supplies inside the BSC
During the operation	Working too close to the front edge/front grill	- Potential hazards to workers, other people (public) and environment	- Work at center area of BSC
	Move arms frequently when removing or introducing objects into BSC	- Disturbance of inflow air velocity	- Avoid rapidly moving your arms in and out of the BSC - Working carefully
	Hazardous biological agents spills	- Contaminate - Worker infected	- Wear the proper PPE
	Use a piece of equipment that creates air turbulence inside the BSC Ex. vortex, centrifuge	- Disturbance of inflow air velocity	- Place an equipment in the back 1/3 of the BSC - Minimise the amount of materials in BSC

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Task	Hazard	Risk	Control
	Handling biohazardous waste	- Contaminate	- Prepare biohazard containers in BSC
After the operation	Inappropriate trash	- Potential hazards to workers, other people (public) and environment	- To dispose of contaminated solid waste, deposit non-sharps into regular biohazard containers and sharps into sharps containers only
	Not decontaminating materials before withdrawing them from the cabinet	- Potential hazards to workers, other people (public) and environment	- Disinfect the materials with 70% alcohol before taking them out
	Forget to remove PPE before leaving work area	- Contaminate	- Clean and disinfect used gloves with 70% alcohol before removing into biohazard containers - Remove PPE before leaving the work area
		- Contaminate	- Allow 5 minutes of operation to purge the system before turn off the BSC

Safety precautions

1. Do not use biosafety cabinets when handling toxic, volatile, or flammable materials.
2. Do not use a gas flame whenever possible as it interferes with airflow.
3. The UV light must be turned off when the cabinet is in use.

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