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Equipment Safety Instruction: Centrifuge

Task	Hazard	Risk	Control
Switch on equipment	Electricity	- Workers electrocute	- Use a properly grounded electrical outlet with the proper voltage and current handling capacity
		- Spill liquid cause an electric shock	- Place Hotplate/Stirrer on flat and level surface - Wear and use PPE
		- Fire and explosion	- Do not store volatile flammable materials near a centrifuge
During the operation	Electricity	- Fire and explosion	- Do not centrifuge flammable or explosive materials - Do not leave a pH Meter unattended when in use
	Unbalanced centrifuge	- Liquid spillage - Damage to the centrifuge or personal injury	- Wear and use PPE - Never overfill centrifuge tubes (do not exceed 3/4 full) - Preparation of a spill kit, first aid kit
	Work with chemical and biological materials	- Injury through exposure to aerosols with chemical and biological hazards	- Wear and use PPE - For biohazardous materials, disinfect outside of tubes prior to removal from biosafety cabinet and loading into rotors/safety cups - When centrifuging infectious materials, wait 10 minutes after the rotor comes to a complete stop before opening the lid - Open sealable tubes/safety cups/rotors inside fume hood or biosafety cabinet, depending on hazard. - Place a biohazard label on the centrifuge

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Task	Hazard	Risk	Control
	Centrifugation at too high a speed for long time	- Thermal, fire and explosion	- Open air conditioning and centrifugal compressors
	-Chipped or cracked glassware	- Injury through exposure to aerosols with chemical and biological hazards	- Ensure that the centrifuge is balanced - Inspect all glassware to be used in the centrifuge (chipped/cracked) - Do not open the lid while the rotor is moving
	The centrifuge is not closed completely	- Damage to centrifuge rotor	- Before centrifuging tubes, ensure that these have been correctly closed in accordance to the instructions
Maintenance		- Contaminate	- Preparing a cleaning schedule - Never clean rotors and associated parts with abrasive wire brushes - Decontaminate the outside of the cups/buckets and rotors before and after centrifugation - Unplugging when not in use

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